


OSHM CAST SPLINT

One-Step Hybrid Mesh



USAGE TREND CHANGE



OSHM CASTING

HYBRID MESH



PROPERTY CHANGES INTO INNOV. CASTING



ONE-STEP SOLUTION [ANYTIME & ANYWHERE]

ALL IN ONE KIT Accelerated and safe Fracture care, require little water, the applications are Easy, Fast and Clean. **Best for On-Site Emergency Response.** **No Electric Saw Required for Removal.**

MAX. COST-BENEFIT



SOLID PROFITABILITY [ADULT/LONG ARM]


ONE-STEP CARE ↓ ECONOMICAL & SIMPLE APPLICATION RAPID TECHNIQUE ADAPTATION ↓ MAX. DAILY CARE OF PATIENTS ↓ NO REMOVAL FOR X-RAYS MAXIMIZE PRODUCTIVITY	PRODUCTS	AVE. WEIGHT	APPLICATION	CARES DURING 2.5 HRS
	OSHM CAST / SPLINT	200 g	3 ~ 10 min	15 - 18 PATIENTS
	PLASTER	4 Kg	45 ~ 55 min	3 - 4 PATIENTS
	FIBER-GLASS	600 g	15 ~ 30 mins	6 - 8 PATIENTS



SAFE AND SOLIDITY, WITHOUT COMPETITION

CLEAR X-RAYS PERMEABILITY and **ECOLOGICAL** post-use treatment
Elasticity before hardening → Safe and simple **Accelerated Total Contact**
Light and **Durable** Harmless to the skin and Best for children patients
3D Mesh Structure → Absorbs impacts and prevents greenhouse effect
 → **Best Rehabilitation**

OSHM CAST of TUBULAR DESIGN → **Ventilation** and **Water Permeability** → **Daily Hygiene** **OSHM SPLINT** with EXCEPTIONAL SOFTNESS → Instantly create anatomical splints → Prevention of Inner Folds



PATENTS & CERTIFICATES



ISO 13485 CERTIFIES COMPANY

- Korean Patent
- Korean Utility model
- USA Patent
- USA FDA Registration

- Register in CE
- International Patent
- Japanese Patent
- Chinese Patent

BEST COST-BENEFIT

- RESOURCES SAVINGS
- REDUCE MEDICAL COMPLICATIONS
- REDUCE TREATMENT TIME
- REDUCE ERROR MARGIN
- EASY WASTE DISPOSAL

(MEDICAL INSTITUTION)

WITHOUT COMPETITION

MAX. BENEFITS

- DAILY CONVENIENCE
- VENTILATED, COMFORTABLE, LIGHTWEIGHT & DURABLE
- WELFARE SATISFACTION
- SAFE FOR DIABETICS & CHILDRENS
- BEST REHABILITATION

(PATIENT)

STEADY USE
IN EMERGENCY

MAXIMUM PROFITABILITY

- INVENTORY ROTATION
- STEADY USE & DEMAND
- IMPROVE PRODUCTIVITY VALUE
- OPTIMIZATION OF PROFITABILITY & OPPORTUNITY COST

(INSTITUTION & SUPPLIER)

MAXIMIZE PRODUCTIVITY

- FAST RESPONSE
- SIMPLE APPLICATION IN 3~10 MIN.
- ONE-STEP CARE (ALL IN ONE KIT)
- ONE-STEP SOLUTION (ALL INCLUDED)
- TOTAL CONTACT
- NO ELECTRICAL SAW REQUIRED

(HEALTH PERSONNEL)

SPECIFICATION OSHM CAST & OSHM SPLINT

TECHNIQUE ADAPTATION IN 3 ~ 5 APPLICATIONS



ITEM	MODEL	SIZES	APPLICATION	
OSHM CAST	1	OSHMC-0225	2" x 25CM	SHORT ARM
	2	OSHMC-0245	2" x 45CM	SHORT OR LONG ARM SHORT OR LONG LEG
	3	OSHMC-0335	3" x 35CM	SHORT OR LONG ARM SHORT LEG
	4	OSHMC-0350	3" x 50CM	SHORT OR LONG ARM SHORT OR LONG LEG
PRE - CUT	5	OSHMC-0370	3" x 70CM	LONG ARM SHORT OR LONG LEG
	6	OSHMC-0435	4" x 35CM	SHORT ARM
	7	OSHMC-0455	4" x 55CM	SHORT LEG
	8	OSHMC-0490	4" x 90CM	LONG LEG
	9	OSHMC-0555	5" x 55CM	SHORT LEG
	10	OSHMC-0590	5" x 90CM	LONG LEG
OSHM SPLINT PRE - CUT	11	OSHMS-0330	3" x 30CM	SHORT ARM / BOOT
	12	OSHMS-0388	3" x 88CM	LONG ARM LONG LEG
	13	OSHMS-0438	4" x 38CM	SHORT ARM / BOOT
	14	OSHMS-0476	4" x 76CM	LONG ARM SHORT OR LONG LEG
	15	OSHMS-0576	5" x 76CM	LONG ARM SHORT OR LONG LEG
	16	OSHMS-5114	5" x 114CM	LONG ARM LONG LEG
OSHM SPLINT ROLL	17	OSHMSR-3450	3" x 450CM	SHORT OR LONG ARM / BOOT
	18	OSHMSR-4450	4" x 450CM	SHORT OR LONG LEG
	19	OSHMSR-5450	5" x 450CM	SHORT OR LONG LEG



MAIN OFFICE

Olympic-Ro 32-Gil, Songpa-Gu, Seoul, Republic of Korea 05547

TEL: +82 70 7521 6116

E-MAIL: generalcrm@mcrops.com

E-MAIL: kwkim@mcrops.com

M.CROPS AMERICA

Hamburgo 9-404, Col. Juárez, Mexico City, Mexico CP. 06600

TEL: +52 55 5703 0292

MOBILE: +52 155 2891 8745